

Current tools:

**Sensobase®****Active search for Web-Thermometers and Web-IO Analog-In devices for sending the measurement values to a database server**

Date	Changes and comments	Revisions
July 25, 2019	<ul style="list-style-type: none"> <li>Problem fixed: Under certain circumstances data stopped being recorded from individual devices.</li> <li>Problem fixed: The filler could sometimes at the end of the reverse search record a solid block of NULL values.</li> <li>Filler can now process logger files with a time base of less than 60s.</li> </ul>	<a href="#">4.45</a> <a href="#">ZIP download</a> (approx. 3.6 MB)
August 29, 2017	<ul style="list-style-type: none"> <li>Restructuring of the inventory list</li> <li>HTTPS support</li> <li>Fixed: previously the Web-IO Digital could only record the first 8 counter states.</li> <li>Smallest zoom range in the diagram is now one hour (previously: one day).</li> <li>Info dialog revised, now also shows version (and status) of the background service.</li> <li>Internal changes in configuring the ODBC data source now also works with ODBC drivers (such as MySQL 3.51.30) which do not correctly implement the SAVEFILE keyword</li> <li>Sensobase can process the password entered in the configuration dialog, thereby eliminating separate entry of login data</li> </ul>	<a href="#">4.40</a> <a href="#">ZIP download</a> (approx. 3.4 MB)
June 14, 2016	<ul style="list-style-type: none"> <li>Chart can be zoomed in as desired</li> <li>Fixed: with very long time periods the X-axis label no longer worked</li> <li>Can now also record the counters on the Web-IO Digital in addition to the sensor values of the Web -IO Analog</li> <li>Better reply in sensor dialog for "identify installed sensors" action</li> <li>Combobox now also suggests 900 sec. (15 min.)</li> </ul>	<a href="#">4.30</a> <a href="#">ZIP download</a> (approx. 2 MB)
October 10, 2011	<ul style="list-style-type: none"> <li><b>New feature</b> for closing gaps after the face in database recording using the data logger in the device.</li> <li><b>New feature</b>: simple chart preview</li> <li>Inventorizing revised: also detects unconnected sensors</li> <li>Export of summary tables now also works with SQLite</li> <li>Combobox now also suggests 900 sec. (15 min.)</li> </ul>	<a href="#">4,11</a> <a href="#">ZIP download</a> (approx. 935 kB)
January 28, 2009	<ul style="list-style-type: none"> <li>Sensobase® now functions as a Windows service</li> </ul>	<a href="#">3.01</a> <a href="#">ZIP download</a> (approx. 460 kB)
February 21, 2008	<ul style="list-style-type: none"> <li>Inventorizing supports longer unit strings (up to 32 characters, like in the database, instead of only 4 as previous)</li> <li>Incorrect description of the database format in the manual corrected ("Devices" -&gt; "Sensors" table, "DeviceID" -&gt; "SensorID" column)</li> </ul>	<a href="#">2.02</a> <a href="#">ZIP download</a> (approx. 216 kB)
03.09.2007	<ul style="list-style-type: none"> <li>Contains a simple SQL console</li> <li>New: Data import from Version 1.xx</li> <li>Change in the database format</li> </ul> <p><b>Note! Old databases are then no longer directly readable (only after data import)</b></p>	<a href="#">2.01</a> <a href="#">ZIP download</a> (approx. 311 kB)
May 2, 2006	<ul style="list-style-type: none"> <li>Current devices added</li> </ul>	<a href="#">1.07</a> <a href="#">ZIP download</a> (approx. 138 kB)
23.09.2005	<ul style="list-style-type: none"> <li>First version of Sensobase®</li> </ul>	<a href="#">1.03</a> <a href="#">ZIP download</a> (approx. 138 kB)

